


Beth Israel
Deaconess
Medical
Center



Harvard
Medical
School

Compressed Sensing in Cardiac MR

Mehmet Akçakaya, Ph.D.

Department of Medicine (Cardiovascular Division)

Beth Israel Deaconess Medical Center

Harvard Medical School

Boston, MA

Motivation

long acquisition time impact :

- spatial and temporal resolution
- spatial coverage
- SNR, CNR
- artifact level (motion, contrast media change)

long CMR exams:

- cardiac MR exams are generally long
- comprehensive exams are becoming more common
- emphasis on patient throughput to reduce cost

Outline

- Brief overview of acceleration methods
- A new CS reconstruction algorithm for high-resolution CMR
- CS-accelerated image acquisition for
 - coronary
 - late gadolinium enhancement
 - cine
 - perfusion
 - other CMR applications

Methods for Acceleration

- Partial Fourier
- Parallel Imaging
- Non-Cartesian Trajectories
- Compressed Sensing
- Spatio-temporal methods
(for dynamic imaging)

Multiple Receiver Coils

- In clinical MRI, receiver coil arrays are used.
- They modulate the intensity of the signal based on their spatial locations.



$$\begin{pmatrix} \mathbf{S}_1 \\ \mathbf{S}_2 \\ \vdots \\ \mathbf{S}_{N_C} \end{pmatrix} = \begin{pmatrix} \mathbf{F}_\Omega \mathbf{C}_1 \mathbf{m} \\ \mathbf{F}_\Omega \mathbf{C}_2 \mathbf{m} \\ \vdots \\ \mathbf{F}_\Omega \mathbf{C}_{N_C} \mathbf{m} \end{pmatrix}$$

\mathbf{C}_k : Coil sensitivity profiles


Parallel Imaging

- Utilizes redundancy in acquisition due to phased-array coils.
- SENSE
 - Image-domain least squares solution
- SENSE and GRAPPA are linear reconstruction methods.
- GRAPPA
 - k-space interpolation
 - Interpolation kernels estimated from center of k-space


Compressed Sensing (CS)

Recently proposed MRI acceleration technique¹
Images are compressible in transform domains



Image Domain



Transform Domain



Wavelet



Finite
Differences


¹ Block et al, MRM, 2007; Lustig et al, MRM, 2007.

CS Reconstruction

F_{Ω} : incoherent k-space undersampling operator

S : measured k-space data (undersampled)

- CS reconstruction



Φ : spatial

Ψ : trajectory


τ : weight of sparsity term

CS in Cardiac MR

- Potential for higher acceleration rates
- Limited use in high spatial-resolution cardiac MR applications
 - Blurring and residual artifacts

Aim

- To use patient-specific and anatomic-specific information for improved reconstruction



- We provide a self-supervised learning framework for 3D cardiac MRI reconstruction

LOST





- Coronary images contain 2D patches of similar signal content at various spatial locations

LOST



- Coronary images contain 2D patches of similar signal content at various spatial locations


LOST



Block Matching

- Construct similarity clusters by block matching¹

matched similar
blocks
(similarity cluster
for the reference
block / voxel)
reference
block
check similarity
with
other blocks




- Clusters are produced for every voxel of the image.

¹ Dabov et al, IEEE TIP, 2007!

Thresholding

- Remove aliasing by using the FFT-sparsity of similarity clusters
 - 1) Hard thresholding in FFT domain
 - Captures l_0 norm of similarity clusters
 - If FFT coefficient $< \tau_{ht}$, set to zero. Otherwise unchanged

- Efficiently remove
- Blurring artifacts




Thresholding

2) Wiener Filtering in FFT domain

- Captures weighted l_2 norm of similarity clusters
- Weight each FFT coefficient, f_k by



$$|f_k|^2 / (|f_k|^2 + T_{\text{wie}}^2)$$

- Reduces blurring artifacts
- Useful in later i weights are reliable




Combining Blocks

- Each block may be in multiple clusters



- Combine these by weighted averaging
 - Intuition: Smaller weight to more aliased blocks

LOST




Implementation

Stage 1: Generate low-resolution estimate from center of k-space

- 1) Adaptively identify similarity clusters ($N_b = 8$)
- 2) Threshold via hard-thresholding

Stage 2: From estimate of stage 1

- 1) Adaptively identify similarity clusters ($N_b = 4$)
- 2) Alternate between hard-thresholding and Wiener filtering




Methods

- Performance evaluation of LOST
 - Targeted coronary MRI
 - Retrospective undersampling
 - Prospective undersampling
- Clinical application of LOST
 - Contrast-enhanced coronary MRI
 - Late Gadolinium enhancement

Retrospective Undersampling

- 1.5T Philips Achieva magnet with 5-channel cardiac coil.
- Right and left coronaries ($N_H = 10$)




- Retrospective undersampling:
 - 50x5 lines in the center
 - **randomly discarding the edges**

¹ van den Berg et al, SIAM JSC, 2008; ² Yang et al, IEEE JSIP, 2010.

LAD/LCX Results

- NAV-gated
- ECG-triggered
- T₂-Prep SSFP
- 1×1×3 mm³ resolution
- Retrospective random under-sampling
- 5-channel cardiac coil



RCA: Right Coronary Artery, LAD: Left Anterior Descending Artery, LCX: Left Circumflex Artery


Performance Evaluation

- Our evaluation shows improved
 - subjective image scores
 - sharpness of the RCA
 - mean square errorwith respect to l_1 minimization and TV regularization
- Next step: implement and evaluate LOST in prospective acquisition


Prospective Undersampling

- Random k-space undersampling in SSFP sequences artifacts

- Radial p




Prospective Undersampling



RCA Results


- NAV-gated
- ECG-triggered
- T₂-Prep SSFP
- 1×1×3 mm³ resolution
- **Prospective random under-sampling** with radial profile order
- 5-channel cardiac coil



LOST allows images to be acquired with 4× acceleration even with 5-channel coil
RCA: Right Coronary Artery, AO: Aortic Root, LV: Left Ventricle, RV: Right Ventricle

LAD/LCX Results

- NAV-gated
- ECG-triggered
- T₂-Prep SSFP
- 1×1×3 mm³ resolution
- **Prospective random under-sampling** with radial profile order
- 5-channel cardiac coil



RCA: Right Coronary Artery, LAD: Left Anterior Descending Artery, LCX: Left Circumflex Artery

Contrast Enhanced (CE) Coronary MRI


- **Advantages¹:**
 - higher SNR
 - higher CNR
- **Disadvantages:**
 - long acquisition time (~10-12 minutes)
 - artifacts due to varying inversion time / contrast washout
- **10 healthy Rapid acquisition is needed**
 - 4-fold acceleration

¹ Bi et al, MRM, 2007; Hu et al, MRM, 2010.


CE WH Coronary MRI Results

- IR-SSFP
- NAV-gated, ECG-Triggered
- bolus 0.2 mmol/kg Gd-BOPTA
- imaging after 2 min of contrast
- resolution = $1.3 \times 1.3 \times 1.3$ mm³
- 4-fold k_y - k_z acceleration
- 5-channel phased-array coil
- 2:50 minutes at 70 bpm,
100% gating efficiency

zero-filled




LOST




CE WH Coronary MRI Results

- IR-SSFP
- bolus 0.2 mmol/kg Gd-BOPTA
- resolution = $1.3 \times 1.3 \times 1.3$ mm³
- 4-fold k_y - k_z acceleration
- 5-channel phased-array coil

zero-filled




LOST





CE WH Coronary MRI Results

zero-filled



LOST



RCA: Right Coronary Artery, LAD: Left Anterior Descending Artery, LCX: Left Circumflex Artery

CE WH Coronary MRI Results

- Subjective quality assessment by two blinded readers in consensus
 - 1 = poor, 4 = excellent ¹
 - Overall score: 3.6 ± 0.5

¹ Kim et al, NEJM, 2001

Late Gadolinium Enhancement

- Late Gadolinium Enhancement (LGE) is used for viability studies
 - Quantification of scar volume and border zone morphology¹
- Higher resolution allows²
 - Identification of small areas of scar
 - Improved gray zone characterization

2D with BH
2x2 mm², 10 mm gap



Higher resolution and shorter
acquisition are desirable


¹ Kim et al., Circulation, 2000, Peter et al., NMR, 2007

3D WH LGE Results

- 46 year old male
- hypertrophic cardiomyopathy
- LGE in myocardium

- IR-GRE
- axial acquisition with WH coverage
- resolution = $1.0 \times 1.0 \times 1.5 \text{ mm}^3$
- 3-fold k_y - k_z acceleration
- 5-channel phased-array coil
- 7 minutes total acquisition time

Low Resolution
 $1.5 \times 1.5 \times 4.0 \text{ mm}^3$





LOST





3D WH LGE Results

Low Resolution
 $1.5 \times 1.5 \times 4.0 \text{ mm}^3$




LOST, Rate = 3
 $1.0 \times 1.0 \times 1.5 \text{ mm}^3$




Motion Correction in LGE¹

- 51 year old female
- pericarditis
- LGE in pericardium

Low-Resolution
 $1.7 \times 1.7 \times 5.0 \text{ mm}^3$



with LOST



- IR-GRE
- resolution = $1.3 \times 1.7 \times 1.7 \text{ mm}^3$
- 5-channel phased-array coil


¹ Moghari et al, MRM, in press

Images courtesy of S. Hong and M. H. Moghari

Cine CMR


- Image the cardiac morphology throughout different phases of the cardiac cycle.
- Used for functional assessment.
 - Quantification of ejection fraction, end-diastolic volume, end-systolic volume, stroke volume, etc
- Accelerated imaging allows
 - Higher spatial or temporal resolution
 - Less breath-holds

2D with BH
2×2 mm², 30 ms temp. res.



Retrospective Cine

Fully-Sampled




Zerofilled (retrospective undersampling, $R = 6$)




Transform Domain LOST

l_1 minimization in x-f space , $R = 6$




reconstruction in x-f space, with transform-domain LOST , $R = 6$




Transform Domain LOST

Fully-Sampled



reconstruction in x-f space, with
transform-domain LOST , $R = 6$



kt-SPIRiT for Retrospective Cine

fully-sampled



kt-GRAPPA ($R=7$)



kt-SPIRiT ($R=7$)




Lai, ISMRM 2010

Courtesy of Drs. P. Lai (GE), M. Lustig (UC Berkeley)

kt-SPIRiT for Prospective Gating

fully-sampled



kt-GRAPPA (R=6)




kt-SPIRiT(R=6)




kt-Group Sparse for Cine

kt group sparse

x - f space



signal support



Training data: 12.5%

fully-sampled



$R = 3$



$R = 5$




$R = 7$




Perfusion CMR

- Images wash-in of contrast media with the blood during the initial pass through myocardium
- Used for assessment of perfusion defects
 - Diagnosis of coronary artery disease
- Accelerated imaging allows
 - Higher spatial or temporal resolution
 - Better coverage

Multislice-2D CMR Perfusion



- Multislice TurboFLASH sequence with 8-fold k_y -t acceleration
- Spatial resolution = $1.6 \times 1.6 \times 8$ mm³
- Temporal resolution = 60 ms
- 3T Siemens Tim Trio, 12-element coil




Courtesy of R. Otazo, D. Sodickson, NYU

3D First-Pass CMR Perfusion

- TurboFLASH sequence
- FOV = $340 \times 340 \times 100$ mm³
- 40 dynamics
- temporal resolution = 220.8 ms
- spatial resolution = $2.7 \times 2.7 \times 8$ mm³
- 16-fold k_y - k_z -t acceleration

- 3 Tesla Siemens Verio
- 32-channel body array




Courtesy of Drs. R. Otazo, D. Sodickson, NYU

Other CMR Applications




Non-Contrast PV-MRA


Fully Sampled




$R = 2$




$R = 4$



$R = 6$



MR Angiography with L₁-SPIRiT




Top: 12.5 fold accelerated blood-pool contrast-enhanced extremity MR angiogram in a 4 year old with a parallel imaging alone (ARC) using a dedicated 32 channel pediatric coil. 750 x 750 x 800 μm resolution. Bottom: L1-SPIRiT reconstruction recovers substantial detail, with quality rivaling a catheter angiogram. The fast scan avoids venous contamination.

Courtesy of Drs. M. Lustig and S. Vasanawalla, Stanford

Real-Time Cine MRI


- True-FISP
- 8-fold acceleration
- Temporal resolution = 42.5 ms
- Spatial resolution = 2.3 mm




Feng et al. ISMRM 2010; 3602

CS in Black Blood CMR

Fully sampled




R = 5



CS

R = 5



auto-calibrated
SparseSENSE

Prieto, ISMRM 2010

Courtesy of Drs. C. Prieto, T. Schaeffter, R. Botnar KCL

Conclusions

- Introduced the LOST algorithm, which uses patient-specific and anatomy-specific information for improved reconstruction
- Provided an overview of applications of CS in CMR

The Future

- More clinical validation needed
- Not available by vendors
- Faster and robust reconstruction needed

§ Potential to significantly accelerate and/or improve CMR image acquisition

Acknowledgements

BIDMC Cardiac MR Center

- Reza Nezafat, Ph. D.
- Warren J. Manning, M.D.
- Michael L. Chuang, M.D.
- Neil M. Rofsky, M.D.
- Vahid Tarokh, Ph.D.
- Tamer A. Basha, Ph.D.
- Peng Hu, Ph.D.
- Mehdi H. Moghari, Ph.D.
- Seunghoon Nam, M. S.

Slide Contribution:

- Michael Lustig, Ph.D.
- Rene Botnar, Ph.D.
- Daniel Sodickson, Ph.D.
- Ricardo Otazo, Ph.D.
- Tobias Schaeffter, Ph.D.
- Muhammad Usman

